

# Onze infuuspompen



Toedieningswijze	Cadd Legacy One	Cadd Legacy PCA	Cadd Legacy Plus	Cadd Solis Vip	Asena GW	Bodyguard 323	Bodyguard 323 Color Vision	Asena GS	Crono S-Pid 50 V2	Crono S-Pid 100	T34L PCA	T34L standaard
Subcutaan		✓		✓	✓	✓	✓	✓	✓	✓	✓	✓
Perifeer	✓	✓	✓	✓	✓	✓	✓	✓			✓	✓
CVC	✓		✓	✓	✓	✓	✓					
PAC	✓		✓	✓	✓	✓	✓					
PICC	✓		✓	✓	✓	✓	✓					
Epiduraal				✓				✓			✓	✓
Spinaal				✓							✓	✓

Doseringsfrequentie	Cadd Legacy One	Cadd Legacy PCA	Cadd Legacy Plus	Cadd Solis Vip	Asena GW	Bodyguard 323	Bodyguard 323 Color Vision	Asena GS	Crono S-Pid 50 V2	Crono S-Pid 100	T34L PCA	T34L standaard
Bolus		✓		✓	✓	✓	✓				✓	✓
Continu ml/uur		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Continu ml/24 uur	✓											
Intermitterend			✓	✓								
PCA		✓		✓							✓	
Stappenprogramma				✓								
Min. loopsnelheid	0,1ml/u	0,1ml/u	0,1ml/u	0,1ml/u	1ml/u	0,1ml/u	0,1ml/u	0,1ml/u	0,1ml/u	0,2ml/u	0,1ml/u	0,1ml/u
Max. loopsnelheid	125ml/u	50ml/u	125ml/u	500ml/u	999ml/u	500ml/u	500ml/u	999ml/u	100ml/u	300ml/u	650ml/u	650ml/u

Veel gebruikt voor	Cadd Legacy One	Cadd Legacy PCA	Cadd Legacy Plus	Cadd Solis Vip	Asena GW	Bodyguard 323	Bodyguard 323 Color Vision	Asena GS	Crono S-Pid 50 V2	Crono S-Pid 100	T34L PCA	T34L standaard
Antibioticum	✓	✓*	✓	✓	✓	✓	✓					
Bloedtransfusie					✓							
Corticosteroid			✓	✓	✓	✓	✓					
Cytostaticum	✓		✓	✓	✓	✓	✓					
Diureticum			✓	✓	✓	✓	✓					
Dobutamine			✓		✓	✓	✓					
Enzymtherapie	✓			✓	✓	✓	✓					
Immunoglobuline				✓	✓	✓	✓		✓	✓		
Pijnmedicatie		✓		✓	✓	✓	✓	✓	✓	✓	✓	✓
TPV				✓	✓	✓	✓					
Vocht	✓		✓	✓	✓	✓	✓					

\* Kan alleen 50ml per uur inlopen

Overige kenmerken	Cadd Legacy One	Cadd Legacy PCA	Cadd Legacy Plus	Cadd Solis Vip	Asena GW	Bodyguard 323	Bodyguard 323 Color Vision	Asena GS	Crono S-Pid 50 V2	Crono S-Pid 100	T34L PCA	T34L standaard
Pompgewicht (g)	392	383	392	595	1.490	365	398	2.500	140	140	400	400
Inclusief tas	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓
Verwisselbare batterij(en)	✓	✓	✓	✓		✓	✓			✓		
Verwisselbare accu				✓		✓	✓					
Niet-verwisselbare accu											✓	✓
Netroom	✓	✓	✓	✓	✓	✓	✓	✓			✓	✓

Alarmpuncties	Cadd Legacy One	Cadd Legacy PCA	Cadd Legacy Plus	Cadd Solis Vip	Asena GW	Bodyguard 323	Bodyguard 323 Color Vision	Asena GS	Crono S-Pid 50 V2	Crono S-Pid 100	T34L PCA	T34L standaard
Batterij/accu leeg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Infusie gereed/pomp gestopt	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Lucht in de lijn	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Occlusiedruk	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reservoir leeg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Systeemfout	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

Doseringshoeveelheid	Cadd Legacy One	Cadd Legacy PCA	Cadd Legacy Plus	Cadd Solis Vip	Asena GW	Bodyguard 323	Bodyguard 323 Color Vision	Asena GS	Crono S-Pid 50 V2	Crono S-Pid 100	T34L PCA	T34L standaard
50 ml cassette	✓	✓	✓	✓								
100 ml cassette	✓	✓	✓	✓								
250 ml cassette	✓	✓	✓	✓								
Infuszak	✓	✓	✓	✓	✓	✓	✓					
10 ml spuit								✓			✓	✓
20 ml spuit								✓			✓	✓
30 ml spuit								✓			✓	✓
50 ml spuit								✓	✓		✓	✓
100 ml spuit										✓		
Fles	✓	✓	✓	✓	✓	✓	✓					